



REF **Sea Simulation**

CE Class I Medical Device

User manual

Distribution mode

Available for direct download at
<http://virtualisvr.com/espace-client/>
Use under licence

 **VIRTUALIS**

Avenue de l'Europe - 34830 CLAPIERS - Tel. 09 80 80 92 91

DESCRIPTION

SEA SIMULATION is an immersive 3D simulation software based on virtual reality technology, meaning a person can be immersed in a digitally created artificial world. The software is used for desensitization in cases of naupathia (sea sickness).

INDICATIONS

Treatment of naupathia (kinetosis), landing sickness.

CONTRAINDICATIONS

Epileptic patients, children under 15 years of age, pregnant women

FOR USE BY

Healthcare professionals: Physiotherapists; Occupational therapists; Neuropsychologists; ENT doctors; Neurologists; PMR doctors (physical medicine and rehabilitation), etc.

Research Centres: CNRS, CHU, INSERM, etc.

WARNINGS AND CAUTIONS

During sessions, stay close to the patient in order to anticipate any loss of balance or discomfort caused by the use of virtual reality.

Define a working area of about 3m² to allow for risk-free movements.

Take a 10 to 15 minute break every 30 minutes of use.

Potential adverse effects are those due to the use of Virtual Reality, namely vomiting, malaise, dizziness, syncope.

The accessories required to use the software may emit radio waves that can interfere with the operation of nearby electronic devices. If you have a pacemaker or other implanted medical device, do not use the product until you have taken advice from your doctor or the manufacturer of your medical device.



Any serious incident should be notified in writing to qualite@virtualisvr.com



Table of Contents

- 1. GENERAL..... 4
 - 1.1. General information on travel sickness..... 4
 - 1.2. Advice for use 8
 - 1.3. Hardware and minimum configuration requirements 9
- 2. USE OF PATIENT MANAGEMENT..... 10
- 3. SEA SIMULATION 12
 - 3.1. Start interface..... 12
 - 3.2. Module field of application 13
 - 3.3. Installing the patient 13
 - 3.4. Session settings 14
 - 3.5. Shortcuts 18
 - 3.6. Data processing 19

1. GENERAL

1.1. General information on travel sickness

General information on the treatment of motion sickness:

The objective here is to provide some basics on motion sickness for those who would like to learn more about it.

These are simple recommendations for use, of course, everyone is free to use this software using their own protocols, to adapt them to their patients, etc... These are ideas, guidelines, that everyone can enrich by sharing them on the Virtualis website Forum.

Re-education principle: Habituation

Objective: to gradually reproduce the triggering conditions

- Therefore, for each case, the triggering conditions must be identified in order to best apply our software on motion sickness (Sea Simulation, Waves (Sway Referenced), Reading (Sway Referenced) Car Simulation).

Trigger mechanism theories:

2 mechanisms:

- Unusual stimulation of peripheral sensors (Habituation without Sensory Conflicts)
- Sensory Conflicts (case of our software: Sea Simulation, Waves (Sway Referenced), Reading (Sway Referenced), Car Simulation).

Unusual stimulation of peripheral sensors depending on:

- **Type of movement:** roll, pitch etc.
- **Stimulation frequency:** movement amplitude

Protocol proposal:

At least reproduce fighter pilot training by using "Coriolis" accelerations, (greatly inspired by: "***The Italian Air Force rehabilitation program for air sickness***", M. Lucertini, V. Lugli, *Acta Otorhinolaryngol Ital* 24, 181-187, 2004).

- Rotations in the chair at different speeds (from 40 to 120°/s), eyes closed, adding head movements in the 3 planes (Coriolis effect).
- We measure the time withstood for each speed in the Clockwise and Anticlockwise directions, and for each movement (Flexion/Extension, Tilts, Rotations)
- The objective is to withstand 2 min. of clockwise and anticlockwise before moving to a higher speed

For example, you can start with slow flexion extension (less than 1 mvt/sec), eyes closed, combined with rotations in the chair, at a constant speed of 40°/s in the clockwise direction. Patients, for example, will withstand one minute and 40 seconds. You can wait a few minutes and, depending on their condition, start again in the anticlockwise direction, where, for example, patients will withstand 48s.




TIP: To try to maintain a constant speed in the chair without an instrument, convert the angular speed into "how long for ½ a rotation" and readjust, stopwatch in hand, every ½ turn

40°/s = 9 sec per turn (360/40=9) => 4.5s per ½ turn, so passages at 4.5, 9, 13.5, 18 etc.

60°/s = 6 sec per turn (360/60=6) => 3s per ½ turn, so passages at 3,6,9,12, etc.

	40 °/S	60 °/S	90 °/S	120 °/C
Rotations				
Right / Left				
Clockwise				
Anticlockwise				
Tilts				
Clockwise				
Anticlockwise				
Flex° / Ext°				
Clockwise	1 min40			
Anticlockwise	48s			

Sensory Conflicts



There are at least 2 categories of conflicts depending on the systems involved and 2 types of conflict for each category.

2 Main Conflict Categories

Depending on the sensory systems involved:

Visual / Vestibular

- the most often described in the literature

SC Canals / Otolithic System

- Vision is not involved (Blind people, eyes closed, may be sensitive)

In each category, there are 2 types of conflict

- **Type 1:**

The 2 sensory systems continuously send contradictory or uncorrelated information

For example: Head movements with visual distortion (binoculars, prism), watching waves from a boat.

- **Type 2:**

One system sends information while the signal expected by the other system is absent

For example: Optokinetics, Imax Cinema, Simulator Sickness (fixed base), etc.

Most highly triggering environments involve more than one type of conflict:

TYPE	CATEGORY	
	VISUAL (A)/VESTIBULAR (B)	CANAL (A) / OTOLITH (B)
1 : A and B signals simultaneously give contradictory information	Watching waves from a ship Use of binoculars in a moving vehicle - Making head movements when vision is distorted by an optical device Reading hand-held material in a moving vehicle	Making head movements while rotating about another axis (Coriolis) / Making head movements in an abnormal acceleration environment, (eg, hyper- or hypogravity) Vestibular disorders (eg, Ménière's disease, acute labyrinthitis, trauma)
2a: A signals are received but expected B signals are absent	OptoKiNetics , Simulator sickness (fixed base) / VIMS / VR Sickness 3D Cinema/Imax sickness/ Circular linearvection	Space sickness Pressure (alternobaric) vertigo Caloric stimulation of semicircular Canals
2b: B signals are received but expected A signals are absent	Looking inside a moving vehicle without external visual reference (cabin of a boat) Reading in a moving vehicle	Low-frequency (< 0.5 Hz) translational oscillation OVAR

Summary table for protocols:

	CATEGORY	
TYPE	VISUAL (A)/VESTIBULAR (B)	CANAL (A) / OTOLITH (B)
1 : A and B signals simultaneously give contradictory information	Waves (sway referenced)	Coriolis Protocol
2a: A signals are received but expected B signals are absent	Sea Simulation Car Simulation	Non-physiological condition
2b: B signals are received but expected A signals are absent	Reading (sway referenced)	Non-physiological condition

1.2. Advice for use

Immersion in Virtual Reality is a powerful tool, especially for stimuli that can induce sensory conflicts

These stimulations have the potential to cause certain disorders: Vasovagal syncope, epileptic seizures, migraines, etc. (Despite a test phase on more than 2000 patients, similarly to previous generation optokinetics, caution is required).

This type of re-education must be undertaken progressively, especially in Virtual Reality where the stimulation is "powerful".

The contraindications are identical: Mainly epilepsy and migraines.

As postural reactions can be spectacular, it is VERY STRONGLY advised to place patients in a safe environment and to stay close to them throughout the session.

It is also recommended to increase the duration and intensity of the stimulation very gradually, after an initial short session to make sure of patients' tolerance to this type of stimulation.

The treatment of motion sickness is by "habituation", so symptoms felt during transport must be re-created very gradually. It is absolutely essential to interrupt the session when the first symptoms

appear, usually "sweating". Accepting that some motivated patients wish to go further would be counter-productive. It is up to the healthcare professional to "dose" the immersion to avoid causing neurovegetative symptoms. This type of symptom may intensify during the hour following the session.


Virtualis declines any liability for any disorders suffered by patients during or after use of its software.

1.3. Hardware and minimum configuration requirements

Hardware required to use the system:

- VR Ready PC
- VR System: HTC VIVE, HTC VIVE Pro or compatible system
- Lighthouse bases (HTC VIVE tracking)
- Thrustmaster T150 Wheel and pedals
- XBOX 360 Controllers
- Dynamic posturography platform (MotionVR) (optional)
- USB HUB

In order to install and use our virtual reality applications, we recommend a configuration equal to or higher than the system requirements:




2. USE OF PATIENT MANAGEMENT

Once connected to the Patient Management software, you arrive on the home page. It is from this home page that you will be able to start your VR software as well as the other Patient Management functions.

The software can be grouped according to criteria such as "Assessment" or "Re-education" and then by pathology type: Neurology, Balance, Functional or Kinetosis.


CINETO VR software contains the following modules: Lift, Reading (sway referenced), Car Simulation, Sea Simulation and Waves (sway referenced).

You can start or switch from one software to another from the home page by clicking the corresponding "Start" or "Protocols" button.




Some software can be started either in **manual mode**, by directly clicking the "Start" button, or in **protocol mode** by clicking the "Protocols" button.

The **manual mode** allows users to choose the settings for each environment. The **protocol mode** offers several sessions with different difficulty levels to test and gradually accustom patients to the VR environment.




Software that is not part of your subscription package is greyed out. If you want to use it, please contact our sales department.



3. SEA SIMULATION

3.1. Start interface




When starting the software in **manual mode** ("Start" button), the opening is made in a start interface, consisting of a module selection menu on the left, a set up area on the right, and an action area at the bottom right.

Depending on the module selected in the left menu, the set up area shows the various possible settings/information.

The general Patient Management menu can be accessed from the start interface by simply clicking the "Back" button located in the action area, or by pressing the "escape" key on the

keyboard.

The module is launched by simply clicking the "Start" button in the action area.



Once this button is pressed, the module starts by taking into account the specified settings.

The Start/Quit buttons allow the environment to be played back or stopped entirely to adapt the experience to the patient's sensations.

Once an environment has been selected, it launches in the headset, and you can see and track what is happening in your patients' headset from the software window.

3.2. Module field of application

This simulation is designed specifically for the desensitization of Naupathia (seasickness).

Most often, through immersion, it makes it possible to physically feel abdominal sensations felt in live conditions on a boat in the waves.

Wave frequency and amplitude (height) are configurable.

The boat is fully steerable using the joystick. Possibility of facing the oncoming waves (less triggering) or from the side, to move forward, turn, reverse etc.

The aim here is to reproduce the "Optokinetic" type Visual-Vestibular conflict, i.e. when Vision perceives a movement and the vestibular apparatus does not. Added to that are "physical" sensations very close to those actually felt on a boat.

The realistic environment adds a psychological immersion aspect in the triggering environment.

Visual-Vestibular Conflict

- Visual Information: Present
- Vestibular Information: Absent (in seated position)

3.3. Installing the patient

To work specifically on seasickness, it is recommended to seat the patient.

However, it is still possible to use the software in a standing position, which causes very strong postural reactions that can lead to falls.

Stay in contact with the patient who is at risk of falling +++ during vertical stimulations, or install them in a closed and safe environment, such as a dynamic posturography platform fitted with a guardrail or safety harness.

This type of work in the standing position can be preferred when stimulating patients' proprioception, sense of balance, etc.

3.4. Session settings

The variable settings for this module are as follows:

Wave frequency

Corresponds to the number of waves per second; the higher this value, the more waves will be generated

Wave height

The higher this value, the higher the waves will be

MotionVR settings

Platform amplitude values can be defined either by using the cursor or by choosing one of the proposed values by simply clicking on the button corresponding to the value.

Display mode

Two possibilities: Headset or Screen mode



You can choose between "headset" or "screen" mode by simply clicking on the corresponding icon; the selected icon turns blue.

Headset mode : The environment launches into the patient's headset

Screen mode : This mode requires having a second screen connected to the computer



Warning: in this mode the headset must not be used, and the patient must be placed in front of the second screen


StaticVR settings

Raw data sent by the platforms

Yellow dots: Center of Pressure (COP) of each foot

Blue dot: Global Center of Pressure (COP)

The weight distribution for each foot is displayed



Smoothed data & settings:

Tare

Platform reset (must be carried out when empty)

Smoothing

Smoothing force applied to the data

Sensitivity

Multiplier applied to received data

Decrease to reduce motion sensitivity

MotionVR settings


Movement settings

Platform amplitude values can be defined either by using the cursor or by choosing one of the proposed values by simply clicking on the button corresponding to the value.

There is the possibility to choose an overall or per-axis movement amplitude, the presets provide a smooth transition.


For example:

- To work in anteroposterior mode, reduce the left and right amplitude
- To work in the mid-lateral position, reduce the forward and backward amplitude




Movement preview

Used to view the platform tilt, height and amplitude settings (the action area is delimited by a black circle).



Options



Platform positions

Park position: forces the platform to ground level

Regular position: allows the platform to move normally, as provided for in the software


Centered position: forces the platform into a horizontal position at its operating height

Activate manual control

Use the arrow keys on the keyboard and the "+" and "-" keys on the numeric keypad (up and down) to move around.

Enable sway referenced mode (for StaticVR platforms)

The platform movement is controlled by the patient's center of gravity

The  button on the launch interface at the bottom right provides access to other advanced options:

Check the corresponding box to "**Enable development options**"

Two types of display are possible:

- [StaticVR] Display graphics
- [StaticVR] Display port status

These boxes are development options used to control devices connected to the computer. It is not advisable to use them (slows down the software).

Recording Options:

You have the possibility of choosing the type of data to be recorded by ticking the corresponding box:

StaticVR:

- Raw data
- Smoothed data

MotionVR

- Gross position (pitch & height)


Select the location on the computer to save the data and press the "Start Recording" button



If you press the emergency stop button, the platform will freeze in its current position.



A window opens and the following message appears



Warning:




Stopping the software or changing the window may cause the platform to reset to the default position and cause movement that could be dangerous for the patient. It is therefore advisable not to touch the computer again until you have secured the patient when the emergency stop is triggered.

3.5. Shortcuts

Keyboard or joystick shortcuts can be accessed in two ways:

- on the "Shortcuts" tab available at the start interface level
- within the module, by clicking on the joystick icon in the upper right corner of the screen



3.6. Data processing

Data retrieval and analysis uses the Patient Management software.